



SRI INDU INSTITUTE OF ENGINEERING AND TECHNOLOGY

Accredited by NAAC A+ Grade, Recognized under 2(f) of UGC Act 1956.

(Approved by AICTE, New Delhi and Affiliated to JNTUH, Hyderabad)

KhalsaIbrahimpatnam, Sheriguda(V), Ibrahimpatnam(M), Ranga Reddy Dist., Telangana – 501510

DEPARTMENT OF HUMANITIES AND SCIENCES

ENGINEERING CHEMISTRY LAB MASTER MANUAL

(R18-IB.Tech)

CERTIFICATE

*ThisistoCertifyThatPreparedbyDr.Premalatha,G.Vijay
V.Swapna,K.Mounika,K.Sharanya,have preparedmanualof
practical work done in the Engineering Chemistry laboratory ,.*

LabIncharge

H.O.D

INDEX

	Experiments	
1	Determination of total hardness of water by complexometric method using EDTA.	13
2	Determination of chloride content of water by Argentometry.	15
3	Estimation of an HCl by Conductometric titrations.	17
4	Estimation of Acetic acid by Conductometric titrations.	20
5	Estimation of HCl by Potentiometric titrations.	23
6	Estimation of Fe^{2+} by Potentiometry using KMnO_4 .	26
7	Determination of rate constant of acid catalyzed hydrolysis of methyl acetate.	29
8	Synthesis of Aspirin and Paracetamol.	33
9	Thin layer chromatography calculation of R_f values. e.g. ortho and para nitro phenols.	35
10	Determination of acid value of coconut oil.	38
11	Verification of Freundlich adsorption isotherm-adsorption of acetic acid on charcoal.	39
12	Determination of viscosity of castor oil and groundnut oil by using Ostwald's viscometer.	41
13	Determination of partition coefficient of acetic acid between n-butanol and water.	44
14	Determination of surface tension of a given liquid using Stalagmometer.	46
	ADDITIONAL EXPERIMENTS	48
15	Estimation of Mohr's salt using standard KMnO_4 solution.	49
16	Determination of salt concentration by Ion-exchange resin.	51

JNTUH (R18) B.Tech-I Syllabus
CH106BS/CH206ES:ENGINEERINGCHEMISTRYLAB

B.Tech.I YearI Sem

LTP C
0031.5

ListofExperiments:

1. Determinationoftotalhardnessofwaterbycomplexometricmethodusing EDTA
2. Determinationof chloride contentofwaterbyArgentometry
3. EstimationofanHClbyConductometrictitrations
4. EstimationofAcetic acidbyConductometric titrations
5. EstimationofHClbyPotentiometrictitrations
6. EstimationofFe²⁺byPotentiometryusingKMnO₄
7. Determinationofrateconstantofacidcatalysedhydrolysis ofmethylacetate
8. SynthesisofAspirinandParacetamol
9. Thin layer chromatography calculation of R_f values. eg ortho and para nitrophenols
10. Determinationofacidvalueofcoconutoil
11. VerificationofFreundlichadsorptionisotherm-adsorptionofaceticacidon charcoal
12. Determination of viscosity ofcastor oil and ground nut oil by using Ostwald's viscometer.
13. Determinationofpartitioncoefficientofaceticacidbetweenn-butanoland water.
14. DeterminationofsurfacetensionofagiveliquidusingStalagmometer.

References

1. Seniorpracticalphysicalchemistry,B.D.Khosla, A.GulatiandV.Garg (R. Chand & Co., Delhi)
2. Anintroduction to practicalchemistry,K.K.SharmaandD.S.Sharma (Vikas publishing, N. Delhi)
3. Vogel'stextbookofpracticalorganicchemistry5thedition
4. TextbookonExperimentsandcalculationsin Engineeringchemistry–S.S. Dara



SRI INDI INSTITUTE OF ENGINEERING & TECHNOLOGY

ENGINEERING CHEMISTRY LAB

LIST OF EXPERIMENTS:

1. Determination of total hardness of water by complexometric method using EDTA.
2. Determination of chloride content of water by Argentometry.
3. Estimation of an HCl by Conductometric titrations.
4. Estimation of Acetic acid by Conductometric titrations.
5. Estimation of HCl by Potentiometric titrations.
6. Estimation of Fe^{2+} by Potentiometry using KMnO_4 .
7. Determination of rate constant of acid catalysed hydrolysis of Methyl acetate.
8. Synthesis of Aspirin and Paracetamol.
9. Thin layer chromatography calculation of R_f values. eg ortho and para nitrophenols.
10. Determination of acid value of coconut oil.
11. Verification of Freundlich adsorption isotherm-adsorption of acetic acid on charcoal.
12. Determination of viscosity of castor oil and ground nut oil by using Ostwald's viscometer.
13. Determination of partition coefficient of acetic acid between n-Butanol and water.
14. Determination of surface tension of a given liquid using Stalagmometer.

ADDITIONAL EXPERIMENTS:

15. Estimation of Mohr's salt using standard KMnO_4 solution.
16. Determination of salt concentration by Ion-exchange resin.

References

1. Vogel's textbook of practical organic chemistry 5th edition
2. Textbook on Experiments and calculations in Engineering chemistry – S.S. Dara

COURSE OBJECTIVES:

1. To bring adaptability to the concepts of chemistry and to acquire the required skills to become a perfect engineer.
2. To impart the basic knowledge of atomic, molecular and electronic modifications which makes the student to understand the technology based on them.
3. To emphasize on the types of hardness of water, disadvantages of using hard water in industries and necessity to soften hard water by various modes.
4. To acquire the knowledge of electrochemistry, corrosion and water treatment which are essential for the engineers and in industry.
5. To acquire the skills pertaining to spectroscopy and to apply them for medical and other fields.
6. To impart the knowledge of stereochemistry and synthetic aspects useful for understanding reaction pathways.

COURSE OUTCOMES:

- C116.1. Student is able to acquire the scientific attitude by means of distinguishing, analyzing and solving various engineering problems.
- C116.2. Student is able to gain the knowledge of atomic, molecular and electronic changes, and band theory related to conductivity.
- C116.3. Student is able to differentiate between hard & soft water, their effects when used in Thermal power plants and various methods involved in water softening.
- C116.4. Student is able to apply the principles and concepts of Electrochemistry, corrosion and mechanism associated with corrosion control methods to produce advanced corrosion resistant materials like PVC, PS.
- C116.5. Student is able to apply the concept of basic spectroscopy to medical and other fields.
- C116.6. Student is able to acquire the knowledge of configurational and conformational analysis of molecules and reaction mechanisms.

InstructionsforRecordWriting

1. Writeontherighthandpagethefollowinginorder:
 - a) Serialnumberanddateofperformance(inthemargin)
 - b) Nameandnumberoftheexperimentasgiveninthelist
 - c) Aimofthe experiment
 - d) Descriptionoftheapparatus
 - e) Procedureincludingourcesoferrorandprecautionstakentoeliminateorto minimize them.
 - f) InferenceorResult.
 - g) Explanation,ifnecessaryofanydivergenceintheexpectedresult.

2. Lefthandpageshouldcontainthefollowingintheirproperplaces.
 - a) Neatdiagramofthemainapparatusused
 - b) Observationintabularform
 - c) Calculationintabularform
 - d) Graphsheetsandotherpaperstobeattached.

- N.B.(a)Studentsshouldsubmitarecordofthepreviousperimentswhenthey comefor practical work.
- (b) Anexperimentisdeemedtobecompletewhenitissatisfactorilyperformed and recorded.
 - (c) Allsuggestionsgivenordefectsnotedshouldbeattendedtowhentherecords are taken back.

KEEPTHERECORDBOOKNEATITFETCHESMARKS

EXPERIMENTS

Experiment-1**Determination of Total Hardness of Water by Complexometric Method Using EDTA**

Objective: To estimate the total hardness of water by EDTA Method.


Outcome: The students will be able to analyse the nature of salts causing hardness and to solve the engineering problems arising during steam production in Boilers.

Chemicals Required: Ammonia buffer solution, Eriochrome Black-T (EBT), Eriochrome Black - T (EBT), Magnesium chloride, EDTA Solution, Hard water


Apparatus: Conical Flask, Pipette, Burette, Beaker

Principle:

Ethylene diamine tetra acetic acid (EDTA) forms stable complexes with Ca^{+2} and Mg^{+2} ions present in water at pH 9-10. The sample of hard water must be treated with ammonia buffer solution and EBT - indicator which forms unstable, wine red coloured complex with Ca^{+2} and Mg^{+2} present in water.



The stability of a metal indicator complex is less than that of metal EDTA complex. During the titration of the complex with EDTA, EDTA extracts the metal ions from the metal - ion - EBT indicator complex and forms stable, colorless complexes by releasing the free indicator. Where by the end point of the titration is the colour change from wine red to blue.



Wine red coloured complex Colour	Colorless & stable complex	Blue in
-------------------------------------	----------------------------	---------

Procedure:**Preparation of standard hard water:**

Dissolve one gram of pure, dry CaCO_3 in minimum quantity of diluted HCl and evaporate the solution to dryness on water bath. Dissolve the residue in small amount of water and transfer into a 100 ml standard flask. Make up the solution to the mark with distilled water and shake the flask well for uniform concentration.

Standardization of EDTA Solution:

Pipette out 10 ml of standard hard water solution into a conical flask, add 2 ml of buffer solution and 2-3 drops of EBT indicator and titrate the wine red coloured complex with EDTA solution taken in the burette, after rinsing the burette with EDTA solution, till the wine red colour of the solution changes to blue colour. Note the burette reading repeat the titration to get concurrent values.

Estimation of hardness of sample water:

Pipette out 2 ml of the water sample into a 25 ml conical flask, add 2 ml of buffer solution and 3 drops of EBT - indicator. Titrate the wine red coloured solution with EDTA taken in the burette, till a clear blue coloured end point is obtained. Repeat the titration to get concurrent values.

Observation & Calculations

Standardization of EDTA:

S.No.	Volume of hard water	Burette reading		Volume of EDTA
		Initial	Final	
1				
2				
3				

Estimation of Hardness:

S.No.	Volume of Sample hard water	Burette reading		Volume of EDTA
		Initial	Final	
1				
2				
3				

Molarity of standard Hard water CaCO_3

$$M_1 = \frac{\text{Wt of CaCO}_3}{\text{M.wt of CaCO}_3} \times \frac{1}{\text{Volume of solution}}$$

$$= \text{_____ M}$$

Molarity of EDTA Solution (M_2) = ?

$$\frac{V_1 M_1}{n_1} = \frac{V_2 M_2}{n_2}$$

$$n_1 = n_2 = 1$$

$$M_2 = \frac{V_1 M_1}{V_2} = \text{_____ M}$$

$$V_1 = \text{Volume of standard hard water} = \text{_____ ml}$$

$$V_2 = \text{Volume of EDTA} = \text{_____ ml}$$

$$M_1 = \text{Molarity of standard hard water} = \text{_____ M}$$

$$M_2 = \text{Molarity of EDTA} = \text{_____ M}$$

Molarity of the sample water = _____ (M_3) or (M_{sample})

$$M_2 V_2 = M_3 V_3$$

$$M_3 = \frac{M_2 V_2}{V_3} = \text{_____ M}$$

$$V_2 = \text{volume of EDTA}$$

$$M_2 = \text{Molarity of EDTA}$$

$$V_3 = \text{volume of hard water}$$

$$\text{Total hardness of water} = M_3 \times 100 \times 1000$$

$$= \text{_____ ppm.}$$

Result:

The total hardness of sample of water = _____ ppm.

Experiment-2**Determination of chloride content of water by Argentometry**

Objective: To determine chloride ion concentration in a water sample by Argentometry (Mohr's method).

Outcome: The students shall be able to identify & estimate the amount of chloride ions present in given water sample.

Chemicals required: Silver nitrate (AgNO_3), NaCl, Potassium chromate (K_2CrO_4) indicator.



Apparatus: Beaker, burette, pipette, conical flask, measuring jar, volumetric flask

Principle:

Chloride ion (Cl^-) is one of the major inorganic anions in water and wastewater. Along the sea coast chloride may be present in high concentration because of leakage of salt water into the sewage system. It also may be increased by industrial process. A high chloride content may harm metallic pipes and structures as well as growing plants. The measured chloride ions can be used to know salinity of different water sources. For brackish water (or sea water or industrial brine solution), it is an important parameter and indicates the extent of desalting of apparatus required. Generally the Mohr's method is used to estimate the chloride present in water sample.

The Mohr method uses chromate ions as an indicator in the titration of chloride ions with a silver nitrate standard solution. After all the chloride has been precipitated as white silver chloride, the first excess of titrant results in the formation of a silver chromate precipitate, which signals the end point (1).

The reactions are:

**Procedure:**

Preparation of 5% K_2CrO_4 (indicator): Dissolve 1.0 g of K_2CrO_4 in 20 ml of distilled water.

Preparation of standard AgNO_3 solution: Transfer 9 g of AgNO_3 to a 500 ml volumetric flask and made up to volume with distilled water. The resulting solution was approximately 0.1 M.

Standardization of AgNO_3 : Transfer 0.2500 g portions of NaCl into conical flask and dissolved in about 100 mL of distilled water. In order to adjust the pH of the solutions, add small quantities of NaHCO_3 until effervescence ceased. Add about 2 ml of K_2CrO_4 and titrate the solution to the first permanent appearance of red $\text{Ag}_2\text{Cr}_2\text{O}_4$.

Determination of Cl^- in solid sample: Take 20 ml of unknown water sample in to conical. Add small quantities of NaHCO_3 until effervescence ceased. Add about 2 ml of K_2CrO_4 and titrate the solution to the first permanent appearance of red Ag_2CrO_4 .

Interpretation of data**Standardization of AgNO₃:**

Replicate	Weight of NaCl (gr)	Volume of AgNO ₃ (ml)	Concentration of AgNO ₃ (M)
1	0.25	V ml	M ₁
2	0.25	V ml	M ₂
3	0.25	V ml	M ₃

Calculations for Replicate 1 of standardization:

Molecular mass of NaCl = 58.44 g/mole

$$\text{Mmoles of AgNO}_3 = \frac{0.25 \text{ gr NaCl}}{58.44 \text{ gr/mole}} \times \frac{1000 \text{ m moles NaCl}}{1 \text{ mole NaCl}} = 4.278 \text{ mmoles}$$

$$\text{AgNO}_3 = \frac{4.278}{V \text{ ml AgNO}_3} = M$$

$$\text{Molarity of AgNO}_3 = \frac{M_1 + M_2 + M_3}{3} = M$$

Determination of Chloride in Unknown:

Replicate	Volume of water sample	Volume of AgNO ₃ (ml)	Concentration of Cl ⁻ solution
1	20 ml	V ml	
2	20 ml	V ml	
3	20 ml	V ml	

Atomic mass of Cl⁻ = 35.45 g/mole

$$\text{Mmoles of Cl}^- = M_{\text{AgNO}_3} \times V_{\text{AgNO}_3} = M \times V = x \text{ Mmoles}$$

$$= x \text{ M moles} \times 35.45 = y \text{ mg}$$

Experiment-3**Estimation of an HCl by Conductometric Titrations**

Objective To determine the strength of the strong acid by titration with strong base Conductometrically.


Outcome: The students shall be able to analyze the variation of conductance values of given acid with addition of strong base using conductivity meter.

Chemicals Required: Sodium hydroxide, Hydrochloric acid

Apparatus: Conductivity Bridge, Conductivity cell, Burette, Beakers, Standard flask, pipette, Burette Stand.

Principle:


At first solution contain H^+ and Cl^- ions. Since H^+ ions possess greater mobility it follows that the conductivity is mainly due to H^+ ions. The addition of NaOH is represented by the equation.



As NaOH is added the H^+ ions are removed. The conductivity decreases as Na^+ ions do not possess much mobility. At the neutralization point the solution contains Na^+ ions and Cl^- ions and will have minimum conductance value. If NaOH is further added this will add OH^- ions and so the conductivity increases.

Procedure:

A standard solution of 1M NaOH is prepared. Similarly, 0.1M HCl is prepared. 200 ml of HCl is taken in a 250 ml beaker. Further, the conductivity cell is washed with distilled water and rinsed with acid solution. The cell is kept in acid containing beaker and it is connected to the bridge. The conductivity of the solution is measured by adjusting the reading. NaOH solution is taken into burette and add 1 ml of solution to acid, stirred well and conductance is measured. Each time 1 ml of base is added to acid stirred well and the conductance is measured. For every instance equal numbers of values are taken on either side of the point of maximum. Repeat the procedure of addition of 1 ml NaOH and noting the conductivity of the resulting solution. Take 20-25 readings



Calculations:Volume of unknown acid (V_1) = 200 ml.

S.No	Volume of NaOH (mL)	Observed conductance
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		

FORMULA:

$$N_1 V_1 (\text{HCl}) = N_2 V_2 (\text{NaOH})$$

$$N_1 = \frac{1 \times V_2}{200}$$

$$\text{Strength of solution} = \text{Eq. Weight} \times \text{Normality}$$

$$= 36.5 \times N_1$$

$$= \text{----- g/lit}$$

RESULT:

The Normality of strong acid (HCl) determined by titrating against a strong base (NaOH)
= _____ N

Precautions:

1. The conductivity cell should be handled very carefully as it is very delicate.
2. Stirrings should be done after each addition of the titrant.

Advantages or applications of conductance titration:

1. This method can be used to very dilute solutions.
2. Gives very accurate endpoints with an error of ± 0.5
3. These titrations are very useful in case of coloured solutions which cannot be titrated by ordinary volumetric method because colour change of indicator is not clear.
4. Useful for titrating weak acids against weak bases, which otherwise do not give sharp endpoints.
5. No keen observation is necessary near the endpoints since it is detected graphically.

Experiment-4**Estimation of Acetic Acid by Conductometric Titrations**

Objective: To determine the strength of the weak acid by titration with strong base Conductometrically.



Outcome: The students shall be able to analyze the variation of conductance values of given weak acid with addition of strong base using conductivity meter.

Chemicals Required: Sodium hydroxide, acetic acid

Apparatus: Conductivity Bridge, Conductivity cell, Burette, Beakers, Standard flask, pipette, Burette Stand.

Principle:

Before titration low initial conductance is observed due to low H^+ obtained during dissociation of weak CH_3COOH . During titration we can observe slight decrease of conductance due to consumption of H^+ . During progress of titration we can observe slight increase in conductance due to the presence of CH_3COO^- & Na^+ and nearly constant H^+ due to the buffer action of the produced CH_3COONa and the remaining CH_3COOH . After end point excess $NaOH$ will lead to increase in conductance due to increasing of Na^+ and OH^- .

**Procedure:**

Prepare a standard solution of 1M NaOH. Similarly, prepare 0.1M acetic acid. 200 ml of acetic acid is taken in a 250ml beaker. Further, the conductivity cell is washed with distilled water and rinsed with acid solution. The cell is kept in acid-containing beaker and it is connected to the bridge. The conductivity of the solution is measured by adjusting the reading. NaOH solution is taken into burette and add 1 ml of solution to acid, stirred well and conductance is measured. Each time 1 ml of base is added to acid stirred well and the conductance is measured. For every instance equal numbers of values are taken on either side of the point of maximum. Repeat the procedure of addition of 1 ml NaOH and noting the conductivity of the resulting solution. Take 20-25 readings

CALCULATIONS:Volume of unknown acid (V_1) = 200 ml.

S.No	Volume of NaOH	Observed conductance
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
15		
16		
17		
18		
19		
20		

FORMULA:

$$N_1 V_1 (\text{CH}_3\text{COOH}) = N_2 V_2 (\text{NaOH})$$

$$N_1 = \frac{1 \times V_2}{200}$$

Strength of solution = Eq. Weight X Normality

$$= 36.5 \times N_1$$

$$= \text{----- g/lit}$$

RESULT:

The normality of weak acid (acetic acid) determined by titrating against a strong base (NaOH) = $_N$

Precautions:

1. The conductivity cell should be handled very carefully as it is very delicate.
2. Stirring should be done after each addition of the titrant.

Advantages or applications of conductance titration:

1. This method can be used to very dilute solutions.
2. Gives very accurate endpoints with an error of ± 0.5
3. These titrations are very useful in case of coloured solutions which cannot be titrated by ordinary volumetric method because colour change of indicator is not clear.
4. Useful for titrating weak acids against weak bases, which otherwise do not give sharp endpoints.
5. No keen observation is necessary near the endpoints since it is detected graphically.

Experiment-5**Estimation of HCl by Potentiometry**

Objective: To determine the equivalence point between strong acid and strong base and to determine the normality of HCl by titrating with NaOH using potentiometer.


Outcome : The student shall be able to analyze the variation of EMF values of given acid with addition of strong base using potentiometer.

Chemicals required: 0.1N HCl, 1N NaOH, distilled water.

Apparatus: Potentiometer, standard cell, saturated calomel electrode, platinum electrode, beaker, burette, stirrer etc.

Principle

When a solution of strong acid (HCl) is titrated with the solution of a strong base (NaOH), the change in pH will be reflected in the change in EMF. When a small amount of alkali is added to the acid, a little change in the EMF is produced in the beginning. This change in electrode potential depends on the fraction of hydrogen ions removed. As an equivalence point is reached, the fraction of the hydrogen ions removed by constant volume of standard alkali increases rapidly, thereby causing a rapid change in the EMF. Thus if the EMF of the cell is plotted against the volume of the standard alkali added, a curve is obtained. As the change in EMF is much more rapid near the equivalent point, the exact equivalent point is obtained by differential method where, a graph of $\frac{\Delta E}{\Delta V}$ Vs. volume of alkali added, gives the maximum of the curve which corresponds to equivalence point of the titration.



The cell can be represented as


**Procedure****Estimation of HCl:**

10 mL of 0.1N HCl is pipetted out into 100 mL beaker and saturated with Quinhydrone and the indicator electrode (platinum electrode) is dipped. The indicator electrode and saturated calomel electrode (reference electrode) are connected to the potentiometer. The two half cells are connected by means of a salt bridge. The potentiometer is standardized and used for measuring the emf directly. 1N NaOH is taken in the burette and is added to the HCl solution.

The emf of the acid, taken in the beaker, is measured initially. Rough titration is first carried out by adding 1 mL of NaOH and the emf for each addition is measured. The emf decreases gradually and then shows a sudden decrease in the emf. Enough readings are taken after the sudden decrease in emf. From these rough titrations, the range of end point is determined.

After finding the endpoint range, fair titrations are carried out by repeatedly adding 1 mL is added and emf is measured for each addition and the readings are tabulated.

Graph is plotted with volume of alkali (NaOH) along X-axis and measured emf along Y-axis. A sigmoid curve is obtained and the equivalence point is noted at the end point of intersection. To obtain a sharp endpoint, another graph of volume of alkali (NaOH) along X-axis and $\frac{\Delta E}{\Delta V}$ along Y-axis is plotted. The maxima obtained in the curve gives the accurate equivalence point.

Observations:**Pilot Titration:**

Volume of alkali (NaOH) (mL)	EMF (mV)
0.0	
1.0	
2.0	
3.0	
4.0	
5.0	
6.0	
7.0	
8.0	
9.0	

Final Titration

Volume of alkali (mL)	EMF (mV)	$\frac{\Delta E}{\Delta V}$
0		
1.0		
2.0		
3.0		
4.0		
5.0		
6.0		
7.0		

8.0		
9.0		
10.0		
11.0		
12.0		
13.0		
14.0		
15.0		
16.0		

Calculations:**Strength of HCl**

Volume of acid taken in beaker, $V_1 = 20 \text{ mL}$

Strength of acid, $N_1 = ?$

Volume of NaOH (from the graph), $V_2 = \dots\dots\dots \text{ mL}$

Strength of NaOH, $N_2 = \dots\dots\dots$

According to Volumetric Law

$$V_1 N_1 = V_2 N_2$$

$$N_2 = \frac{V_1 N_1}{V_2}$$

Amount of HCl

Amount of HCl present in the given solution = $\dots\dots\dots (N_2) \times 36.45 \text{ g}$

RESULT:

1. The equivalence point of potentiometric titration between strong acid vs strong base is ... mL
2. The Normality of HCl by titrating with NaOH using potentiometer is $\dots\dots\dots \text{ N}$.
3. The amount of HCl present in the given solution = $\dots\dots\dots \text{ g}$.

Experiment-6**Estimation of Fe^{2+} by Potentiometry using KMnO_4**

Objective: To estimate the Fe^{2+} by Potentiometry using KMnO_4 .

Outcome : The student shall be able to analyze the variation of EMF values of given acid with addition of KMnO_4 using Potentiometer.


Chemicals: KMnO_4 , Distilled water, sulphuric acid, Ferrous sulphate solution.

Apparatus: volumetric flask, ptelectrode, saturated calomelelectrode, potentiometer

Principle:

Potentiometric titration is the titration in which potentiometric measurements are carried out in order to fix the end point. In this method, the interest is with the change in electrode potential, rather than with an accurate value for the electrode potential in a given solution. In a potentiometric titration, the change in cell e.m.f. occurs most rapidly in the neighbourhood of the end point.

The $\text{Fe(II)}-\text{KMnO}_4$ redox system is represented as




The determining factor is the ratio of the concentrations of the oxidised and the reduced forms of the iron species.

For the reaction,



$$E = E^\circ + \frac{0.0591}{n} \log \frac{[Ox]}{[Red]}$$

where E^0 is the standard Reduction Potential of the system. Thus the potential of the immersed electrode is controlled by the ration of these concentrations. During redox reactions, the potential changes more rapidly at the vicinity of the end point. The indicator electrode is usually a bright platinum wire or foil, the oxidising agent is taken in the burette. The cell can be represented as,



Here Pt is the indicator electrode and calomel is the reference electrode.

Procedure:

Preparation of 0.1N KMnO_4 :

0.1 N KMnO_4 is prepared by dissolving 0.31 g of analar crystals in distilled water in a 100 ml volumetric flask. The solution is made up to the mark.


Calibration of The Potentiometer:

A standard cell of known emf is connected to the instrument and its emf is set in the voltage scale. The galvanometer key is pressed to complete the circuit and the deflection of the galvanometer needle is noted. If there is any deflection, the current passing through the rheostat is adjusted for null deflection. This procedure makes sure that the value of emf which is read on the scale is the true potential of the cell considered. The potentiometer is calibrated using the Weston standard cell of potential 1.018 V.

Estimation Of Fe (II):

- The given Fe(II) solution is made up to 100 ml in volumetric flask.
- 20 ml of the solution is pipetted out into a clean beaker. To this, 25 ml of 2.5M H_2SO_4 and 50 ml of distilled water are added.

- A platinum electrode is dipped into this solution, and it is coupled with a calomel electrode through a salt bridge. The resulting cell is connected to the potentiometer.
- Standard KMnO_4 solution is added from the burette, to this solution, in steps of 1 ml and the emf is recorded after each addition.
- At the end point, there is a jump in emf due to the absence of Fe^{2+} . The approximate range of the end point is determined.
- The experiment is repeated by adding the titrant in steps of 0.1 ml near the endpoint. A graph is plotted between emf, E and the volume of dichromate added. The inflexion point gives the volume of titrant at the end point.
- The first derivative ($\Delta E/\Delta V$ vs. Volume of titrant) curves give the exact volume of dichromate required for the reaction. From the plot of E vs. Volume of titrant, potential at the equivalent point is obtained.
- Atomic weight of Fe is 55.85

**Table:**

S.No	Volume of KMnO_4	EMF	ΔE	$\Delta E/\Delta V$	V_{mean}
1					
2					
3					
-					
-					

Volume of $\text{KMnO}_4 =$

V_1 Concentration of $\text{KMnO}_4 = N_1 = 0.1 \text{ N}$

Volume of Fe^{2+} solution = $V_2 = 20 \text{ ml}$

Concentration of Fe^{2+} solution = $N_2 = ?$

According to volumetric law $N_1 V_1 = N_2 V_2$

$$N_2 = N_1 V_1 / V_2$$

$$\begin{aligned} \text{Amount of Fe}^{2+} \text{ ion present in the given solution} &= (N_2 \times \text{Equivalent weight of Fe}^{2+}) / 10 \\ &= (N_2 \times 55.85 / 10) \text{ gm} \end{aligned}$$

Result: The amount of Fe^{2+} ion present in the given solution = gm

Experiment-7

Determination of Rate Constant of Acid Catalyzed Hydrolysis of Methyl Acetate

Objective: To determine the rate constant of the hydrolysis of Methyl acetate using an acid as a catalyst.


Outcome: The students shall be able to identify and analyze the rate constant values of methyl acetate hydrolysis at regular intervals of time.

Chemicals: Methyl acetate, distilled water, HCl, NaOH, Phenolphthalein indicator, 0.1 N oxalic acid.

Apparatus: Conical flask, volumetric flask, burette, pipette, reagent bottles.

Principle:

Methyl acetate undergoes hydrolysis, in the presence of an acid (HCl, for example), to give acetic acid and methyl alcohol.



In the presence of an acid, this reaction should be of second order, since two molecules are reacting. But, it is found to be first order. This may be explained in the following way: The rate of the reaction is given by

$$\frac{dx}{dt} = k^1 [\text{CH}_3\text{COOCH}_3] [\text{H}_2\text{O}]$$

where k^1 is the rate constant (or specific rate constant). Since water is present in large excess, its active mass (molar concentration) virtually remains constant during the course of the reaction. Therefore, its active mass gets included in the constant, and the above equation reduces to:

$$\frac{dx}{dt} = k^1 [\text{CH}_3\text{COOCH}_3]$$

Thus, the rate of the reaction is determined by one concentration term only (that is, by a single power of the concentration term only). Hence, the reaction is first order. Such reactions are also referred to as pseudo first order reactions. The progress of the reaction (hydrolysis of ester) is followed by removing a definite volume of the reaction mixture, at definite intervals of time, cooling it in ice, and titrating the acetic acid formed against alkali, which has already been standardized. The amount of alkali used is equivalent to the total amount of hydrochloric acid present originally and the amount of acetic acid formed in the reaction.

The amount of acetic acid formed (x), at definite intervals of time (t), can be obtained. The amount of acetic acid formed, at the end of the reaction, is equivalent to the initial concentration (a) of the ester. Suppose the volumes of the sodium hydroxide solution (titre value) required for neutralization of 5 ml of the reaction mixture are:

- (i) at the commencement of the reaction is V_0
- (ii) after time (t) is V_t
- (iii) at the end of the reaction is V_∞ . The

n:

x (amount of acetic acid formed after time) is proportional to $(V_t - V_0)$. a

(initial concentration of ester) is proportional to $(V_\infty - V_0)$.

$[a-x]$ (concentration of ester present after time t) is proportional to $(V_{\infty} -$

$$V_0) - (V_t - V_0) = (V_{\infty} - V_t)$$

The first order rate expression given by:

$$k_1 = \frac{2.303}{t} \log \frac{a}{[a-x]} \quad \text{would correspond to :}$$

$$k_1 = \frac{2.303}{t} \log \frac{(V_{\infty} - V_0)}{(V_{\infty} - V_t)}$$

Hence, the rate constant (k_1) could be calculated.

Procedure

Step I:

Standardization of NaOH using standard Oxalic acid (0.1N)

1. Pipette out 10 ml of given 0.1N standard Oxalic acid into a 100 ml conical flask.
2. Titrate this solution against the given unknown concentration of NaOH using phenolphthalein indicator until the end point is colorless to pale pink.
3. Tabulate the values and repeat the titration for concurrent readings and determine the unknown concentration of supplied NaOH solution.

Table 1

S.No	Volume of Oxalic Acid taken (ml)	Burette Readings (ml)		Volume of NaOH consumed (ml)
		Initial	Final	
1				
2				
3				

$$N_1 V_1 = N_2 V_2$$

Here, N_1 = concentration of oxalic acid, V_1 = Volume of oxalic acid N_2 =

concentration of NaOH, V_2 = Volume of NaOH

Step II:

Standardization of HCl using NaOH solution

1. Pipette out 2 ml of given HCl into a 100 ml conical flask.
2. Titrate this solution against the NaOH using phenolphthalein indicator until the endpoint is colorless to pale pink.
3. Tabulate the values and repeat the titration for concurrent readings and determine the unknown concentration of supplied HCl solution.

Table 2

S.No	Volume of HCl taken (ml)	Burette Readings (ml)		Volume of NaOH consumed (ml)
		Initial	Final	
1				
2				
3				

$$(\text{NaOH})N_2V_2 = N_3V_3(\text{HCl})$$

Concentration of HCl, $N_3 =$ _____

Step III:

Determination of rate constant (k_1) for the acid-catalyzed hydrolysis of methyl acetate

1. Take 100 ml of given HCl (whose strength is determined in step II) solution in a stoppered reagent bottle.
2. Add 5 ml of methyl acetate solution to the HCl solution. Note the time when half of the methyl acetate solution is added. The mixture is shaken well.
3. Pipette out 5 ml of the reaction mixture and discharge it into 50 ml of ice cold water kept in a conical flask.
4. Titrate the reaction mixture against NaOH solution using phenolphthalein as indicator. This titre value corresponds to V_0 .
5. Repeat steps 3 and 4 at intervals of 5, 10, 15, 20, 30, 45, 60 minutes. Each titre value corresponds to V_t .
6. Take the remaining solution in a stoppered conical flask and heated to 60°C , and keep at this temperature for 5 minutes.
7. The solution is allowed to cool to room temperature.

Table 3

S.No	Time (min)	Volume of solution taken (ml)	Burette Readings (ml)		Volume of NaOH consumed (ml)	$(V_\infty - V_t)$ (ml)	$\log(V_\infty - V_t)$	$k_1 = \frac{2.303}{t} \log \frac{(V_\infty - V_0)}{(V_\infty - V_t)}$ (min^{-1})
			Initial	Final				
1	0 (V_0)							
2	5 (V_{t1})							
3	10 (V_{t2})							
4	15 (V_{t3})							
5	20 (V_{t4})							
6	30 (V_{t5})							
7	45 (V_{t6})							
8	60 (V_{t7})							
9	V_∞							

8. Repeat Steps 3 and 4. This titre value corresponds to V_{∞} till concurrent values are obtained.
9. Plot a graph of $\log(V_{\infty} - V_t)$ versus time (t) and determine the slope.
10. Report the theoretical and graphical value of rate constant (k_1).

Observations and Calculations:

Room Temperature = °C; V_{∞} = ml; V_o = ml; $(V_{\infty} - V_o)$ = ml
;

$\log(V_{\infty} - V_o)$ =

Mean Value of Rate Constant (k_1) = _____

Results:

1. Strength of NaOH Solution = _____

2. Strength of HCl Solution = _____

3. Rate Constant (k_1) for the acid-catalyzed hydrolysis of methyl acetate at ... °C =
= _____ (theoretical).

= _____ (graphical).

Experiment-8


Synthesis of Aspirin and Paracetamol

Objective: To synthesize some common pain relievers: aspirin and Paracetamol.


Outcome : The student shall be able to synthesize the crude sample of Aspirin and Paracetamol.

Principle: Aspirin (acetylsalicylic acid) is both an organic ester and an organic acid. It is used extensively in medicine as a pain killer (analgesic) and as a fever-reducing drug (antipyretic). When ingested, acetylsalicylic acid remains intact in the acidic stomach, but in the basic medium of the upper intestinal tract, it hydrolyzes forming the salicylate and acetate ions.

Aspirin (molar mass of 180.2 g/mol) is prepared by reacting salicylic acid (138.1 g/mol) with acetic anhydride (molar mass of 102.1 g/mol). Aspirin is a weak monoprotic acid.



Paracetamol (acetaminophen) acts as a fever reducer and pain reliever. Acetaminophen is an amide, a compound that is a derivative of ammonia that has been reacted with an acidic substance, in this case, acetic acid. It can be found in several analgesic preparations, such as Tylenol, some of which may contain other ingredients such as caffeine and buffers.



PROCEDURE

a) Preparation of Aspirin

Chemicals: Salicylic acid, acetic anhydride, sulfuric acid, ethanol, ice.

Apparatus: Dropper, Erlenmeyer flask 125 mL, Beakers, Graduated cylinders, Watch glass, Stirring rod, Ring stand, Buchner funnel, filter paper to fit Buchner funnel, vacuum filtration flask, Rubber tubing for vacuum flask, thermometer, dropper.

Procedure

- Weigh out 2.0 g of salicylic acid. Place it in a 125 mL conical flask. Add 5 mL of acetic anhydride.
- Swirl the flask to wet the salicylic acid crystals. Add 5 drops of concentrated sulfuric acid, to the mixture.
- Gently heat the flask in a boiling water bath for about 10 minutes.
- Remove the flask from the hot water bath and add 10 mL of deionized ice water to decompose any excess acetic anhydride.
- Chill the solution in an ice bath until crystals of aspirin no longer form, stirring occasionally to decompose residual acetic anhydride.

- If oil appears instead of a solid, reheat the flask in the hot water bath until the oil disappears and again cool.
- Setup a vacuum filtration apparatus. Wet the filter paper in the Buchner funnel with 1-2 ml of distilled water.
- Turn on the water aspirator. Decant the liquid onto the filter paper, minimizing any transfer of the solid aspirin.
- Add 15 ml of cold water to the flask, swirl, and chill again. Pour the liquid and the crystals of aspirin onto the filter paper. Repeat until the transfer of the crystals to the vacuum filter is complete.
- Determine the mass of the crude aspirin crystals.

b) Preparation of Paracetamol (Acetaminophene)

Chemicals: Acetic anhydride, Phosphoric acid, Ethanol

Apparatus: Dropper, Erlenmeyer flask, 125 mL Beakers, Graduated cylinders, Watch glass, Stirring rod, Vial to hold aspirin sample, Ring stand Clamp (to hold 125-mL conical flask), Buchner funnel, Filter paper to fit Buchner funnel, Vacuum filtration flask, Rubber tubing, Ice, Dropper.

Procedure:

- Fill a 400 mL beaker about half full with water. Place the beaker and water on a hot plate and bring to a boil.
- Weigh out 1.5 g of p-aminophenol and transfer it into a 125 mL conical flask. (Avoid contact with skin. You may wish to wear gloves.)
- Add 25 mL of water. Add 20 drops of concentrated phosphoric acid (H_3PO_4), and swirl the flask until all of the amine dissolves. If not, add a few more drops of phosphoric acid.
- Turn off the hot plate. Place the flask in the hot water. Carefully add 2 mL of acetic anhydride to the flask. Leave the flask in the warm water for 10 minutes.
- Remove the flask and place it in an ice-water bath. Stir the mixture to crystallize the acetaminophen. You may need to scratch the walls of the flask to start the crystallization. If no crystals appear, add a small seed of acetaminophen to start the crystal formation. Allow the flask to stay in the ice-water bath for 30 minutes.
- Collect the crystals in a Buchner funnel using vacuum filtration. Wash the crystals with 10 mL of cold water. Allow the crystals to dry.
- Determine the mass of the crude acetaminophen.

Results:

Mass of salicylic acid = gms

Mass of crude aspirin = gms

Yield of aspirin = %

Mass of P-aminophenol = gms

Mass of paracetamol = gms

Yield of paracetamol = %

Experiment-9

Thin Layer Chromatography Calculation of R_f Values. Eg Ortho and Para Nitro Phenols

Objective: To calculate the R_f values of ortho and para nitrophenols by using thin layer chromatography.

Outcome: The students shall be able to identify and calculate the R_f values of various compounds present in the given mixture.

Chemicals: O-nitrophenol, P-nitrophenol, ethyl acetate, di-chloromethane

Apparatus: TLC plate, capillary tubes,

Principle:

Most reactions produce more than one product. Naturally occurring materials are only rarely 100% pure. It is therefore desirable to have a simple, fast and efficient way to determine the purity of organic mixtures. The separation of a mixture by passing it, in solution, over an adsorbent (such as Alumina or Silica Gel) is the basic idea of Chromatography. It involves the passage of a mobile phase across a stationary phase in a column. Usually a mixture of compounds is present in the mobile phase. As soon as the mixture comes in contact with the stationary phase, some or all of the components of the mixture are adsorbed on it. As additional mobile phase comes along, some or all of the mixture will dissolve and continue moving. This adsorption/solution process continues along the length of the column. If a proper choice of mobile phase, stationary phase, solvent and other operating parameters was made, the mixture will be separated in the column and its various components will emerge at different times.

In Thin Layer Chromatography, a liquid solution is directly applied to a solid adsorbent. Capillary action draws a developing solvent up the TLC plate. As this solvent passes through the spot, the mixture will be dissolved and will begin to move with the solvent front. However, the adsorbent will also absorb part or all of the mixture. As more solvent comes by, the mixture will again go into solution, move further and be reabsorbed. Since different materials will be dissolved and reabsorbed at different rates, separation will take place. This passage of the solvent front through the adsorbent is known as **developing** the plate. The extent of separation, measured by retention factor ("R_f") value differences, will depend on the relative solubilities and relative strengths of adsorption of the components of the mixture.

Organic compounds interact with adsorbents by a variety of interactions. If the compound is non-polar, it can only have weak 'Van der Waals' attractions for the adsorbent. However, more polar molecules may interact more strongly by a variety of mechanisms including dipole-dipole interactions, coordination, and hydrogen bonding. The most important rule of chromatography is that *the more polar compounds will be absorbed most strongly on adsorbents (stationary phases), while non-polar compounds will be only very weakly absorbed*. In a typical chromatography experiment, the non-polar compounds, since they are poorly absorbed, will be held least strongly and will move quickly through the plate. Polar compounds, on the other hand, will be slowed on their process through the plate by their strong interactions with the solid phase. This separation based on polarity will explain most of the chromatography encountered in this course.

Types of Adsorbents used in Chromatography

Listed in decreasing power of adsorption:

Alumina > Activated Charcoal > Magnesium Silicate > Silica > Starch

Solvents Commonly Used in Chromatography


Listed in decreasing polarity:

Acetic Acid > Water > Methanol > Acetone > Ethyl acetate > Diethyl ether > Chloroform > Methylene chloride > Toluene > Cyclohexane > Petroleum ether

For a typical separation, a variety of different combinations of solvent and adsorbent may be effective. Once you have developed your plate, it must be **visualized**. This visualization may be accomplished by reacting the developed plate with a chemical reagent. Iodine (I_2) is one of the easiest to use of these several common chemical visualizing agents. The developed slide is simply exposed to I_2 vapors in a chamber similar to the developing chamber for a few minutes. Almost all compounds will form a weak colored complex with the I_2 . This complex will appear as a darker area on the slide. The 'spots' are characterized by their R_f value, a measure of how far the spot traveled with that combination of adsorbent and solvent.

Procedure:

- Take 1 TLC plate handle it only on the edges, as fingerprints contain UV-active materials.
- Using a pencil draw a very light line across the sheet (short dimension) about 1 cm from one end.
- Then make 4 small light marks at even intervals along the line for spotting the samples. Draw another light line about 1 cm from another end of the plate for the solvent front.
- Obtain a TLC chamber and place solvent, a 5% ethyl acetate in dichloromethane to 0.5 cm height.
- Place a piece of filter paper around the inside surface of the container and extend into the solvent.
- A glass jar with a lid or a beaker with a watch glass or a cover of a Petri dish can be used as a TLC chamber:
- Using clean capillary tubes carefully spot four samples at two pencil marks as shown below.
- The spots should be as small as possible in order to minimize tailing and overlapping when the TLC plate is developed. If a more intense spot is desired, let the spot dry and re-spot in the same location.
- When the spots are dry, place the TLC plate in the developing chamber. Then gently close the chamber.




A: o-Nitrophenol

B: p-Nitrophenol

- Be sure that the bottom edge of the TLC plate is in the solvent but the spots are above the solvent, and the filter paper does not touch the chromatographic sheet.
- When the solvent has moved to the front line, remove the plate. Lay it on a clean surface in a well ventilated area and allow the solvent to evaporate until the plate appears dry.
- Visualize the plate under iodine chamber. Measure all the distances traveled by the compounds and solvent.

- Calculate the retention factor (R_f) for each compound and immediately draw a light pencil line around each spot.

Calculation of retention factor:



$$R_f = \frac{\text{distance travelled by the component}}{\text{distance travelled by the solvent}}$$

$$= \frac{a}{b}$$

Retention factor $R_f = \frac{\text{Distance travelled by the component}}{\text{Distance travelled by the solvent}}$

Result:

Retention factor of o-nitrophenol =

Retention factor of p-nitrophenol =

Experiment-10

Determination of Acid Value of Coconut Oil

Objective: To Determine the Acid Value of Coconut Oil.

Outcome: The students shall be able to identify and calculate the acid values of coconut oil.

Chemicals: NaOH, distilled water, Phenolphthalein indicator, potassium hydrogen phthalate, coconut oil, ethanol, diethyl ether

Apparatus: Volumetric flask, conical flask, pipette, burette, weighing balance, beaker

Principle:

Acidity may occur in coconut oil upon storage due to the decomposition of the oil to free fatty acids, peroxides and low molecular weight aldehydes and ketones, which produce distinctive smell and affect the quality of oil. Acid value may be defined as the number of milligrams of KOH required to neutralise the free acid in 1 gm of fatty oil.

Procedure:

Preparation of Solutions:

- 1) **0.1N Sodium hydroxide:** Weigh accurately 0.4 gm of sodium hydroxide into 100 mL standard flask, dissolve in little distilled water and make up the solution to the mark with distilled water and shake the flask well for uniform concentration.
- 2) **Ethanol-ether solution:** Prepare a mixture of ethanol and diethyl ether (1:1 v/v); add 1 mL of phenolphthalein indicator and titrate with NaOH till pale pink colour is observed. This is the neutralisation of ethanol-ether solution.
- 3) **Standard solution of potassium hydrogen phthalate:** Accurately weigh about 1.2 gm of potassium hydrogen phthalate into 100 mL flask dissolve in little water make up the solution to the mark with distilled water.

Procedure:

- 1) **Standardisation of NaOH solution:** Pipette out 20 mL of potassium hydrogen phthalate solution into a 250 mL conical flask and add 2 drops of phenolphthalein indicator and titrate the solution with NaOH taken in burette. Pale pink colour is the end point of the titration. Note the burette reading and let the value be x mL
- 2) **Standardisation of coconut oil:** Accurately weigh 1 gm of coconut oil into a 250 mL conical flask; add 50 mL of ethanol-ether solution. Shake the flask well for uniform dissolution of the sample. Add 1 drop of phenolphthalein indicator and titrate the solution with NaOH taken in the burette till pale pink colour is observed. Note the burette reading, let the value be y mL

Calculations:

$$(M_1) \text{ Molarity of potassium hydrogen phthalate} = \frac{1.2 \times 1000}{204.22 \times 100} = 0.058 \text{ M}$$

$$(M_2) \text{ Molarity of NaOH: } M_1 V_1 / M_2 V_2$$

$$M_2 = M_1 V_1 / x$$

$$\text{Acid value of coconut oil} = \frac{M_{\text{NaOH}}}{W} \quad (W = \text{weight of coconut oil} = 1 \text{ gm})$$

Result: The acid value of given oil is

Experiment-11

Verification of Freundlich Adsorption Isotherm-Adsorption of

Acetic Acid on Charcoal

Objective: To study the adsorption of acetic acid from aqueous solution by activated charcoal & examine the validity of Freundlich & Langmuir's isotherm.

Outcome : The student shall be able to identify and calculate the adsorption values & plot the adsorption isotherms.

Chemicals: 0.1N acetic acid solution, 0.1N sodium hydroxide solution, phenolphthalein indicator, powdered activated charcoal etc.

Apparatus: Five reagent bottles with stoppers, burette, pipette, conical flasks

Principle: Adsorption is accumulation of a substance at an interface. The adsorption of a solute from a solution, generally follows the Freundlich empirical adsorption isotherm given by

$$\frac{x}{m} = KC^n$$

Where

x = weight of adsorbent, m = mass of adsorbent, K = constant representing the capacity of the adsorbent, and C = equilibrium concentration of the solution.

From the concentration (C), the volume of each original solution and concentration of each solution after equilibrium (C_e), the weight of oxalic acid per gram of the adsorbent (x/m) is calculated as follows.

$$\frac{x}{m} = \frac{(C - C_e) \times 60}{10}$$

Draw a graph of $\log(x/m)$ vs. $\log C_e$ and calculate the constant n and K using the equation.

$$\log \frac{x}{m} = \log K + \frac{1}{n} \log C_e$$

Procedure:

1. Take five clean, dry stopper reagent bottles & label them from 1 to 5.
2. Add by means of a burette 50, 40, 30, 20 & 10 ml of 0.1N acetic acid solution & 0, 10, 20, 30 & 40 ml distilled water in bottle nos. 1, 2, 3, 4 & 5 resp.
3. Weigh accurately 1 gm of activated charcoal & add to bottle no. 1. Similarly, add 1 gm activated charcoal in each of the remaining bottles. Stopper the bottles & shake them well.
4. Meanwhile determine the initial concentration of acid by titrating against std. NaOH solution using phenolphthalein indicator.
5. Using dry filter papers of dual size filter the solution of each bottle in separated dry flasks. Titrate 10 ml of filtrate from each bottle with 0.1N NaOH using phenolphthalein indicator. Take three readings for each bottle & take the mean. These readings give the equilibrium concentrations of the solution i.e. concentration of acid after adsorption. Tabulate the result.

Observation table:

Bottle	Initial	Equilibrium	Amount of	x/m	$\log x/m$	$\log C_e$

No.	conc. of acid C_0 m eq./lit	Conc. of acid, C_e gm eq./lit	acid adsorbed (x)			
1						
2						
3						
4						
5						

Here $x = \frac{(C_0 - C_e) \times V \times \text{Eq. wt. of the acid}}{1000}$

Where, V = total volume (50ml) of the solution in each bottle.

Calculation:

Calculate initial concentration (C_0) of acetic acid in gm equivalent per lit for each sample. Also calculate the equilibrium concentration (C_e) of the acid solution in each of the bottles in gm equivalent per lit. Then calculate the amount of the acid adsorbed in each bottle as follows:

Graph:

1. Plot the graph of $\log(x/m)$ (Y-axis) against $\log C_e$ (X-axis). A straight line will be obtained. This is in agreement with Freundlich equation. Find slope & Y-intercept values.
- 2) Plot $C_e/(x/m)$ (Y-axis) against C_e (X-axis). A straight line obtained shows the agreement with Langmuir adsorption isotherm. Calculate slope & Y-intercept.

Result:

Thus Freundlich & Langmuir adsorption isotherm are studied by the above experiment.

Experiment-12

Determination of Viscosity of Castor Oil and Ground Nut Oil by

Using Ostwald's Viscometer

Objective: To determine the viscosity of the given castor oil and groundnut oil.

Outcome: The students shall be able to determine and calculate the viscosity coefficient values of Castor oil and groundnut oil.

Chemicals: castor oil and groundnut oil, distilled water.

Apparatus: Specific gravity bottle, viscometer, rubber tube with screw pinch cock, stand, beaker

Theory: The force of friction which one part of the liquid offers to another part of the liquid is called viscosity. For measuring the viscosity coefficient, Ostwald viscometer method is used which is based on Poiseuille's law. According to this law, the rate of flow of liquid through a capillary tube having viscosity coefficient ' η ' can be expressed as

$$\eta = \frac{\pi \cdot r^4 t P}{8 v l}$$

where, v = vol. of liquid (in ml)

t = flow time (in sec.) through capillary r =

radius of the capillary (in cm)

l = length of the capillary (in cm)

P = hydrostatic pressure (in dyne/sq. cm)

η = viscosity coefficient (in poise).

Since, the hydrostatic pressure (the driving force) of the liquid is given by $P = d g h$ (where h is the height of the column and d is the density of the liquid);

$$\eta \propto P t; \text{ or, } \eta \propto d g h t$$

If, η_1 and η_2 are the viscosity coefficients of the liquids under study, d_1, d_2 are their densities and t_1 and t_2 are their times of flow of *equal volume* of liquids through the same capillary respectively, then

$$\eta_1 \propto d_1 g h t_1$$

$$\eta_2 \propto d_2 g h t_2$$

Hence

$$\frac{\eta_1}{\eta_2} = \frac{d_1 t_1}{d_2 t_2}$$

Here, usually the viscosity of given liquid is measured with respect to water whose viscosity is known very accurately at different temperatures. The SI physical unit of viscosity is the pascal-second (**Pa·s**), (i.e., $\text{kg} \cdot \text{m}^{-1} \cdot \text{s}^{-1}$). This means: if a fluid with a viscosity of one **Pa·s** is placed between two plates, and one plate is pushed sideways with a shear stress of one pascal, it moves a distance equal to the thickness of the layer between the plates in one second. The cgs unit for the same is the **poise (P)**, (named after J. L. Marie Poiseuille). It is more commonly expressed, as **centipoise (cP)**. [1 cP = 0.001 Pa·s]. Water at 20 °C has a viscosity of 1.0020 cP.

Procedure:

1. Note the laboratory temperature.
2. Wash the density bottle with distilled water and dry.
3. Take the weight of the empty & filled (with distilled water) specific gravity bottle (with stopper). Then, weigh the filled with specific gravity bottle unknown given liquids

individually. Use the data for measuring the densities.

4. Clean and rinse the viscometer properly with distilled water. Fix the viscometer vertically on the stand and filled with specific amount (say 20ml) of mixture (every time take the same volume).

5. Time of flows was recorded for each solution (water and the given liquids).

6. Take 3 to 4 readings.

Observations:

1. Laboratory temperature = °C

2. Density measurement:

Weight of empty density bottle (w₁) = ... g.

Weight of density bottle with water (w₂) = ... g.

Weight of density bottle with castor oil (w₃) = ... g.

Weight of density bottle with groundnut oil (w₄) = ... g.

So,

Weight of water = (w₂ - w₁) = ... g.

Weight of castor oil = (w₃ - w₁) = ... g.

Weight of groundnut oil = (w₄ - w₁) = ... g.

Sample	Flow times			
	t ₁	t ₂	t ₃	Average
Water				
Castor oil				
Groundnut oil				

Calculations:

Determination of the viscosity of the liquid (η)

As discussed earlier

$$\frac{\eta_c}{\eta_w} = \frac{d_c t_c}{d_w t_w}$$

$$\eta_c = \frac{d_c t_c}{d_w t_w} \times \eta_w$$

here η_c = viscosity of castor oil; η_w = viscosity of water (1.0020 cP at 20°C) d_c =

density of castor oil; d_w = density of water

t_c = flow time for castor oil; t_w = flow time for water

$$\frac{\eta_g}{\eta_w} = \frac{d_g t_g}{d_w t_w}$$

$$\eta_g = \frac{d_g t_g}{d_w t_w} \times \eta_w$$

here η_c = viscosity of groundnut oil; η_w = viscosity of water

d_c = density of ground nut oil; d_w = density of water

t_c =flowtimeforgroundnut; t_w =flowtimeforwater

Result:

Theviscosityof thegivencastoroilwithrespectto wateratlaboratorytemperaturewas found to be..... CP.

Theviscosityof thegivengroundnutoilwithrespect towateratlaboratorytemperature was found to be..... CP.

Experiment-13

**DeterminationofPartitionCoefficientofAceticAcidbetween n-
Butanol and Water**

Objective: To determine the partition coefficient of acetic acid between n-butanol and water.

Outcome: The students shall be able to determine and calculate the Partition coefficient values of n-butanol and water.

Chemicals: Acetic acid, n-butanol, distilled water

Apparatus: Stopper bottle, beaker, pipette, conical flask

Principle:

In dilute solutions at constant temperature a solute which exists in the same molecular species in two non-miscible solvents, will distribute itself between these two solvents at constant temperature according to the partition law, the partition coefficient:

$$k = C_1/C_2$$

where

C_1 = concentration of acetic acid in water, and

C_2 = concentration of acetic acid in butyl alcohol

The concentrations C_1 and C_2 must be expressed in the same unit either as grams, gram molecules, or gram equivalents per litre.

Procedure:

- Boil about 200 mL of distilled water in a beaker for 10 minutes. Pour into a flask and stopper lightly and cool. This is CO_2 free water for later use.
- In a 200 mL glass stoppered bottle place about 70 mL of approximately 2M acetic acid and 50 mL of n-butyl alcohol.
- Stopper the bottle and shake well for at least 1 minute, and then allow the liquid layers to separate. Note and record the temperature of the mixture.
- Insert a 25 mL pipette and carefully withdraw a 25 mL aliquot of the upper alcohol layer. The pipette should first be rinsed by sucking up a little of the solution and discarding this.
- Pipette the solution into a second glass stoppered bottle and add an approximately equal volume of boiled distilled water to this second bottle.
- Shake well to transfer the acid to the water layer, add 3 drops of phenolphthalein and titrate with the approximately 1.0 M sodium hydroxide. The bottle should be stoppered from time to time and vigorously shaken and titration continued until a faint permanent pink colour remains.
- Pipette, also a 25 mL aliquot from the lower aqueous layer of the first bottle as follows. Close the pipette with the finger and place it carefully in the lower layer.
- Suck up and blow out gently a small quantity of liquid to wash out any small quantity of the upper layer that has got into the pipette. Rinse the pipette with a little of solution.
- Allow the solutions to settle and withdraw 25 mL of the lower solution. Place this solution in a flask, add three drops of phenolphthalein and titrate with the approximately 1.0 M sodium hydroxide.
- Add about 25 mL each of fresh butanol and boil distilled water (but no further acid) to the original mixture remaining in the first bottle. Repeat the procedure above and sample and titrate with the new concentrations of acid.

The following table summarises the above procedure:

Initial mixture: 50 mL n-butyl alcohol, 70 mL 2M acid:

Summary of procedure

Alcohol Layer		Water Layer		Alkali for Titration mL		K = C ₁ /C ₂
Sample Removed ml	Fresh Alcohol Added	Sample Removed ml	Water Added ml	Water Alcohol C ₁	C ₂	
25	--	25	--			
--	25	--	25	--	--	
25	--	25	--			

Result: The partition coefficient of acetic acid between n-Butanol and water is.....

Experiment-14

Determination of Surface Tension of a Given Liquid Using Stalagmometer

Objective: Determine the surface tension of a given liquid at room temperature using stalagmometer by drop number method.

Outcome : The student shall be able to determine and calculate the surface tension values of reference liquid and given liquid.

Apparatus: Stalagmometer, specific gravity bottle, a small rubber tube, screw pinch cork

Chemicals: Distilled water, experimental liquid

Principle: In the drop number method, the number of drops formed by equal volumes of two liquids is compared. If m_1 and m_2 is the mass of one drop of each of the liquid having densities d_1 and d_2 respectively. If n_1 and n_2 is the number of drops formed by volume v of the two liquids, then their surface tensions are related as

$$\gamma_1/\gamma_2 = (d_1/d_2) * (n_2/n_1)$$

One of the liquids is water its surface tension and density are known. Then the surface tension of the given liquid can be calculated.

Procedure:

1. Clean the stalagmometer with chromic acid mix, wash with water and dry it
2. Attach a small piece of rubber tube having a screw pinch cock at the upper end of the stalagmometer.
3. Immerse the lower end of the stalagmometer in distilled water and suck the water 1-2cm above mark A. adjust the pinch cork so that 10-15 drops fall per minute .
4. Clamp the stalagmometer allow the water drops to fall and start counting the number of drops when the meniscus crosses the upper mark A and stop counting when the meniscus passes mark B
5. Repeat the exercise to take three to four readings
6. Rinse the stalagmometer with alcohol and dry it
7. Suck the given liquid in the stalagmometer and count the drops as in case of water
8. Take a clean dry weighing bottle weigh it with water as well as with liquid.
9. Note the temp of water taken in a beaker.

Observations:

Room temp = $t^{\circ}\text{C}$

Density of water = d_w

Surface tension of water = γ dynes/cm

No of drops From a Fixed Volume				Mean
Liquid	1....	2.....	3.....	$n_l =$
Water	1....	2.....	3.....	$n_w =$

Weight of empty specific gravity bottle = w_1 gram

Weight of specific gravity bottle + water = w_2 gram

Weight of empty sp. gravity bottle + liquid = w_3 gram

Weight of water = $(w_2 - w_1)$ gram

Weight of liquid = $(w_3 - w_1)$ gram

Calculations:

Densityoftheliquid

$$D_l = (w_3 - w_1) / (w_2 - w_1) * d_w$$

Surfacetensionof liquid=

$$\gamma_l = (d_l / d_w) * (n_w / n_l) * \gamma_w$$

Result

The surfacetensionof liquid is.....dynes/cm.

ADDITIONAL EXPERIMENTS

EXPERIMENTNO-15**ESTIMATIONOFMOHR'SSALT USINGSTANDARDKMnO₄SOLUTION**




Objective: To estimate the amount of Ferrous ammonium sulphate present in one litre solution using 0.02M KMnO₄ solution.

Outcome: The students shall be able to identify the oxidation of Ferrous ion in the given reaction.

Apparatus: Burette, pipette, conical flask, beakers, burette stand, glazed tile, wash bottle.

Chemicals: Mohr's salt solution, 0.02M KMnO₄ solution, dilute sulphuric acid.

Principle: Potassium permanganate in presence of dilute sulphuric acid oxidizes ferrous sulphate to ferric sulphate. Ammonium sulphate does not take part in the reaction.



2mole of KMnO₄ = 10moles of Mohr's salt

Formulae:

$$M_1 V_1 / n_1 = M_2 V_2 / n_2$$

Where M₁ = Molarity of KMnO₄ solution.

V₁ = volume of KMnO₄ solution

n₁ = number of moles of KMnO₄ M₂

= Molarity of Mohr's salt

V₂ = volume of Mohr's salt

n₂ = number of moles of Mohr's salt

Amount of Ferrous ammonium sulphate present in one litre of given solution

= Molarity of Mohr's salt x Gram Mol. Wt of Mohr's

salt

Procedure :

1. Rinse and fill the burette with the given KMnO₄ solution.
2. Pipette out 20ml of Mohr's salt solution into a clean conical flask.
3. Add equal volume of dilute sulphuric acid.
4. Titrate the solution against 0.02M KMnO₄ solution till pale pink colour appears.
5. Note the vol. of KMnO₄ used. It is called End point.
6. Repeat the titration until the concurrent readings are obtained.
7. Calculate the molarity of Mohr's salt solution and amount of Mohr's salt present in one litre of solution by using the above mentioned formulae.

Observations & Calculations:

S.No	Vol. Of Mohr's salt solution	Burette reading		Vol. of KMnO ₄ rundown (b-a) ml
		Initial (a)	Final (b)	
1				
2				
3				
4				

$$M_1 V_1 / n_1 = M_2 V_2 / n_2$$

Where M₁ = Molarity of KMnO₄ solution.

V₁ = volume of KMnO₄ solution

n₁ = number of moles of KMnO₄ M₂ =

Molarity of Mohr's salt

V₂ = volume of Mohr's salt

n₂ = number of moles of Mohr's salt

Amount of Ferrous ammonium sulphate present in one litre of given solution
 = Molarity of Mohr's salt \times Gram Mol. Wt of Mohr's salt
 = ----- \times 392 = ----- g/litre.

Therefore the amount of Mohr's salt present in one litre of given solution = ----- g/litre.

Result :

The amount of Mohr's salt present in one litre of given solution = ----- g/litre.

VIVA QUESTIONS:

1. What is Mohr's salt solution?
2. Write the oxidation reaction involved in the experiment?
3. What is the oxidation state of Chromium ion in $K_2Cr_2O_7$?
4. Define normality?
5. Mention the conditions maintained during titration.

EXPERIMENTNO-16**DeterminationofSaltConcentrationbyIonexchangeResin**

Objective: To determine the NaNO₃ salt concentration by Ion-exchange Resin.

Outcome: The students shall be able to identify the ion exchange reactions & role of ion exchange resins in purification of hard water.

Apparatus: - Separatory funnel-250ml, column chromatographic tubes-2, glass wool or sintered glass disc, beaker, conical flask, Pipette.

Chemicals Required: Cation Exchange Resin, Anion exchange Resin, NaNO₃, NaOH-0.1N, Phenolphthalein, AgNO₃-0.1 N, Potassium Chromate

Principle:

Many substances both neutral and artificial have ion-exchange properties. All ion exchangers have common properties.

- They are almost insoluble in water and organic solvents.
- They contain active or counter ions that will exchange reversibly with other ions in the surrounding solution without any appreciable physical change in the material.
- The ion-exchanger is of complex nature, in fact polymeric that carries an electric charge that is exactly neutralized by the charges on the counterion. If the active ions are cations it is called cation exchanger.

For example widely used cation exchanging resin is obtained by the copolymerization of styrene and divinylbenzene followed by sulphonation to give cation exchanging resin.

The general representation is RH^+ where C^+ = cation.



The anion exchanger is a polymer that contains amine or quaternary ammonium or OH⁻ groups in the interstitial parts of polymer lattice which is generally represented by R^1OH and the exchange takes place as follows

**Procedure:**

Weigh out exactly 0.5gms air dried and transfer into a column with the help of small camel-hair brush, through a funnel into the column. Add sufficient distilled water to cover the resin. Dislodge the air bubbles that stick to resin by applying intermittent pressure to the rubber tubing. Adjust the level of the outer tube, so that the liquid in the column will drain 1 cm above

resin beads. Take a 250 ml separatory funnel with 100 ml of the given NaNO_3 salt solution and run this solution at a rate of 2 ml per minute. Collect the effluent with standard 0.1 N NaOH using phenolphthalein as indicator. Pink color is the end point of the solution.



Take 1 gm anion exchange in a column and take 100 ml of the given sodium nitrate solution into a separatory funnel and add slowly NaNO_3 salt solution at rate of 2 ml per minute. Collect the effluent and titrate with 0.1 N Silver nitrate solution using Potassium chromate as indicator

Result:-

The concentration and amount of NaNO_3 salt = _____

Observation & Calculations:-

Cation exchanger column:-

$$N_1 V_1 = N_2 V_2$$

V_1 = Volume of Std NaOH =

N_1 = Normality of Std NaOH =

N_2 = Normality of Na effluent = V_2

= Volume of Na effluent =

$$N_2 = \frac{N_1 V_1}{V_2}$$

Amount of Na^+ = $N_2 \times 23 =$ _____ gms/100ml

$$N_1 V_1 = N_3 V_3$$

V_1 = Volume of Std AgNO_3 = N_1 =

Normality of Std AgNO_3

= N_2 = Normality of NO_3^- effluent =

V_2 = Volume of NO_3^- effluent =

$$N_3 = N_1 V_1 / V_3$$

Amount of NO_3^- = $N_3 \times 62 =$ _____ gms /100ml

VIVAQUESTIONS:

1. What is Ion-exchanger resin?
2. In estimation, what indicators are employed?
3. What is the principle involved in given experiment?
4. What is the colour change at endpoint of cation exchange effluent?